



MIFILL

Cross-Linked HA Filler


JFDA Approved

Distributed by Balanced Chemicals Trading Co.
Amman-Jordan +962 6582 3429 +962 79 660 7233

• MiFill TECHNOLOGY

PATENTED *MCL

(MULTI-STAGED CROSS LINKING) TECHNOLOGY




Conventional Filler that claimed to be long lasting HA Filler, made by increasing the cross-linking rate. However, due to the limitations of the Cross-Linking technology, most of it only manage to increase the concentration and particle size, without contribution to the longevity of the product. And when a product simply increased on its concentration will give characteristic which is not easily moldable.

MiFill designed with new technology of MCL (Multi-Staged Cross Linking), to modify the existing Cross-Linking technology by several unique steps, created new innovation HA Filler which provides excellent viscosity and moldability.


Cohesiveness & Viscosity

Preventing
“spread out” of the filler
: Migration Resistant




• MiFill DATA ANALYSIS

Low Extruding Force for Smooth injection




Refractive Absorbance Analysis



The wavelength pattern of MiFill is the same as its raw material, showing that the product has not altered its physical properties.

MiFill is safe with no chemical modifications in HA even after its multi-stage cross-linking process.

Hyaluronidase Test : Melting down over-dosed or Incorrectly injected Area



• MiFill APPLICATIONS



HANDS REJUVENATION

MiFill LIGHT

• MIFILL B (BODY FILLER)

PRODUCT DESCRIPTION & INDICATIONS

MiFill B contains biocompatible Cross-Linked Hyaluronic Acid, made from a non-animal derivatives source. It is a pure and stable gel form which will be naturally degraded (biodegradable) in our body.

INDICATIONS

- ✓ Body volume restoration.
- ✓ Body surface contouring & enhancement.
- ✓ Ideal for Breasts enhancement / augmentation.
- ✓ Calves & Buttocks contouring.
- ✓ Smoothens uneven skin surfaces or discrepancies, which is commonly caused by liposuction procedure and any other surgical procedures.

EFFECTS & SAFETY

MiFill B can lasts up to 12-14 months. Since it is made by Hyaluronic Acid which naturally found in our body, MiFill B is biocompatible, biodegradable, and safe. It will be slowly metabolized in our body and gradually disappeared through the normal excretion system mechanism.

Retouch and re-treatment is highly recommended to maintain the ideal long lasting result. The other benefits as the composition made by Hyaluronic Acid, beside giving result for the main purpose of the treatment, it will also stimulates tissue, helps to repair texture, and improves the skin quality.



• MiFill PRODUCTS LINE-UP

Mifill	LIGHT	DEEP	VOLUME	B
Figure				
Compostion	(Cross-Linked) HA 20mg/ml Non Lidocaine	(Cross-Linked) HA 20mg/ml Non Lidocaine	(Cross-Linked) HA 20mg/ml Non Lidocaine	(Cross-Linked) HA 20mg/ml Non Lidocaine
Volume	1.0ml	1.0ml	1.0ml	10.0ml
Proticte Size	200-300um	300-400um	390-480um	400-500um
Duration	6-8 months	10-12 months	12-14 months	12-14 months
Cross Linking rate	●●●○	●●●●	●●●●●	●●●●●
Indication	Forehead Periorbital Lines Under Eyes Lips Neck Lines Hand Rejuvenation	Temples Malar Smiling Lines Cheek Lips Deep Wrinkles	Nose Cheek Chin	Breasts Buttocks Body Filler Labia majora



 **beaumed®**